



# WF User Access Toolkit

LabVIEW addon

This toolkit from WireFlow allows an application to easily grant access to specific code in a LabVIEW application\*, depending on the currently active user. It supports user interactive authentication as well as by security tokens, i.e. login dialog, WF Security dongle (WF 2007 or WF 2008) or the WF Fingerprint reader (WF 2111).


In the configuration tool each user can be configured to use one or more of the authentication methods, allowing an application to mix authentications within one single configuration.

Users can be associated with one or more groups, where each group identifies additional features that is enabled for the group members.

At runtime the application checks if the current user is an administrator, if the user is member of a named group or if a specific feature is enabled.

## Application areas

- User level restrictions
- Feature control
- Headless login



## Features

- Prevent unauthorized access to VIs, menus etc.
- Support for different authentication control methods;
- User-Password
- WF Dongle
- WF Fingerprint reader
- Support headless user authentication
- Support LabVIEW Real-Time targets when headless credentials are used.
- Same interface for all types of authentications
- AES-256 encryption of configuration data

\*requires LabVIEW 2012 or greater

WireFlow AB

Theres Svenssons gata 10  
SE-417 55 Göteborg  
Sweden

[www.wireflow.se](http://www.wireflow.se)

User Access Toolkit  
Data Sheet  
AB0005-040, rev B